



384-well MatriPlate

Dimension	ANSI Standards	Measured
A – Left Edge to A1	11.42 - 12.84	12.13
B – Top Edge to A1	8.28 - 9.70	8.99
C – Length Footprint	127.51 - 128.01	127.76
D – Left Edge to A24	114.92- 116.34	115.63
E – Width Footprint	85.23 - 85.73	85.48
F – Top Edge to P1	75.78 - 77.20	76.49
Plate Height	14.10 - 14.60	14.35
Well Depth (from plate surface to well bottom)		7.83
Range Surface Flatness (+/- 0.075)		_____
Range Well Bottom Flatness (+/- 0.075)		_____

* All measurements in millimeters