


MatriTubes

96 and 384

- Manufactured with high-grade, clarified polypropylene.
- Enables individual compound retrieval and increases compound stability.
- Optional 2D data matrix barcodes.
- Optional MatriCap (96-well)
- MiniTube locking feature (patent pending).
- Both 96 and 384 products offer locking and cassette stacking without the requirement for lids.



96 well format 2D matrix bar codes



96 MiniTube Products

MSP096-T: 5-500uL Bag of natural 96-well individual MatriTubes, Polypropylene

MSP096-2DT: 5-500uL Bag of natural 96-well individual 2D MatriTubes, Polypropylene

MSP096-C: 96-well black MatriTube empty cassette, Polystyrene

MSP096-A2DT: 5-500uL 96 MatriTube Plate filled with 2D barcoded tubes

MSP096-AT: 5-500uL 96 MatriTube Plate filled with tubes

MSP096-MC1: 96-well MatriCap Seal, solid pierceable TPE cap; 12 pierces with 16 gauge probe

MSP096-MC2: 96-well MatriCap Mat, split septum pierceable TPE cap

96 Tube working volume:

Septum seal - 350uL

Heat seal - 450uL

384 MiniTube Products

MSP384-C: 384-well black MatriTube empty cassette, Polystyrene

MSP384-A: 0.5-40uL natural 384 well individual MatriTube storage plate filled with 384 tubes



www.matrical.com